




**CHEM 0012**

**Lab 4: Determination of an Equilibrium  
Constant using Spectroscopy**



- 
- Determination of the equilibrium constant of the following equilibrium system at room temperature.




Iron (III) cation

thiocyanate ion

iron thiocyanate complex


# Experimental Procedure Overview



- Equilibrium concentrations of product and reactant will be determined from five different starting points.
- The equilibrium concentration of the red-brown product will be determined using a spectrophotometer.
- The equilibrium concentrations of the reactants will be calculated.

# Determine Solution Concentrations by Spectrophotometry

## Spectrophotometric Analysis



Measure concentration of a Solute by its absorbance of light.

- Nature of the solute
- Wavelength of light
- Distance light travels through the solution
- Concentration of the solute

where

- $\epsilon$  - molar absorptivity
- $\ell$  - distance light travels through the solution
- $c$  - concentration of the absorbing solute

Plot of "Absorbance versus Concentration" is a straight line of slope  $\epsilon\ell$ . This is called a **Beer's-Law** plot.



# Experimental Procedure - Steps

- Part A

- Prepare 6 solutions to determine the calibration curve for  $[\text{Fe}(\text{SCN})^{2+}]$ .
- Measure the absorbance at wavelength,  $\lambda = 447 \text{ nm}$ .
- Plot a calibration curve and get the equation of the best fit line.

- Part B

- Prepare 5 solutions + 1 blank and allow equilibrium to establish from 5 different starting points.
- Measure the absorbance of the 5 equilibrium solutions at wavelength,  $\lambda = 447 \text{ nm}$
- Determine the equilibrium  $[\text{FeSCN}^{2+}]$  from the calibration curve from Part A.
- Determine the equilibrium  $[\text{Fe}^{3+}]$  and  $[\text{SCN}^-]$  from the ICE table.
- Calculate K

# Prepare Solutions

- Work with partners
- Half the groups start with Part A and half with Part B
- $\text{Fe}(\text{NO}_3)_3$  - source of  $\text{Fe}^{3+}$ 
  - from bottle-top dispenser
  - 0.150 M  $\text{Fe}(\text{NO}_3)_3$  in Part A
  - 0.00150 M  $\text{Fe}(\text{NO}_3)_3$  in Part B
- KSCN – source of  $\text{SCN}^-$ 
  - from bottle-top dispenser
  - 0.000500 M  $\text{SCN}^-$  in Part A
  - 0.00300 M  $\text{SCN}^-$  in Part B
- Prepare 12 solutions – (2 blanks)
  - Total volume 25.00 mL using volumetric flasks
  - Follow recipe for Part A and Part B
  - Top up with 0.10 M nitric acid






Measure Absorbance at wavelength=447 nm

- Solutions transferred in cuvettes. Check cuvettes for:
  - Fingerprints
  - Scratches
  - Air bubbles
  - Use one cuvette only rinsing each time with solution to be analyzed

## Part A - Determine absorbance of standard solution



- 25.00 mL volumetric flask
- Using the 0.150 M  $\text{Fe}(\text{NO}_3)_3$ 
  - $[\text{Fe}^{3+}] > 300 [\text{SCN}^{-}]$
  - Shifts equilibrium almost entirely to right.
  - Essentially all  $\text{SCN}^{-}$  converted to  $\text{Fe}(\text{SCN})^{2+}$
- Absorbance at  $\lambda = 447 \text{ nm}$
- Measure blank with only  $\text{Fe}(\text{NO}_3)_3$  in acid.



# Determine absorbance of all equilibrium solutions

- 10 solutions + 2 blanks
- Absorbance at  $\lambda = 447 \text{ nm}$
- Use one cuvette only rinsing each time with solution to be analyzed
- Check cuvettes for:
  - Fingerprints
  - Air bubbles
  - scratches

## Discard of waste solutions:

- Use a 400 mL beaker and discard all liquid waste at your bench top.
- Discard your waste beaker content in the waste container provided.

# Calculations:

## Part A – Calibration Curve

$$[\text{Fe}(\text{SCN})^{2+}] = \frac{(\text{Volume SCN}^- (\text{L}))(\text{Concentration SCN}^- (\text{M}))}{\text{Total volume (L)}}$$


- Volume  $\text{SCN}^-$  (L) = refer to Table 4-1
- Concentration  $\text{SCN}^-$  (M) = refer to Table 4-1
- Total volume (L) = Volume of volumetric flask in liter

### Calibration Curve:

Plot Absorbance Versus  $[\text{Fe}(\text{SCN})^{2+}]$

- Use the net absorbance column
- Include (0,0) on the graph
- Get the best fit line from Excel

# Calibration Curve



Used to determine  $[\text{FeSCN}^{2+}]$  in equilibrium solutions in Part B

## Calculations:

### Part B - Equilibrium concentration of $\text{Fe}(\text{SCN})^{2+}$

- For Solutions 8 – 12,  $[\text{Fe}(\text{SCN})^{2+}]$  is determined from the calibration curve.  
This is  $[\text{Fe}(\text{SCN})^{2+}]_{\text{equilibrium}} = \text{'x'}$
- Use the Net Absorbance values for each solution.

# Calculations:

## Initial concentration of Reactants

$$\text{Initial [Fe}^{3+}] = \frac{(\text{Volume Fe(NO}_3)_3 \text{ (L)}) (\text{Concentration Fe(NO}_3)_3 \text{ (M)})}{\text{Total volume (L)}}$$

- Volume Fe(NO<sub>3</sub>)<sub>3</sub> (L) = refer to Table 4-2
- Concentration Fe(NO<sub>3</sub>)<sub>3</sub> (M) = refer to Table 4-2
- Total volume (L) = Volume of volumetric flask in liter

$$\text{Initial [SCN}^-] = \frac{(\text{Volume SCN}^- \text{ (L)}) (\text{Concentration SCN}^- \text{ (M)})}{\text{Total volume (L)}}$$

- Volume SCN<sup>-</sup> (L) = refer to Table 4-2
- Concentration SCN<sup>-</sup> (M) = refer to Table 4-2
- Total volume (L) = Volume of volumetric flask in liter

# Calculations:

## Equilibrium concentration of Reactants

	$[\text{Fe}^{3+}]$	$[\text{SCN}^-]$	$[\text{Fe}(\text{SCN})^{2+}]$
[ I ]	$[\text{Fe}^{3+}]_{\text{initial}}$	$[\text{SCN}^-]_{\text{initial}}$	0
[ C ]	- X	- X	+ X
[ E ]	$[\text{Fe}^{3+}]_{\text{initial}} - \text{X}$	$[\text{SCN}^-]_{\text{initial}} - \text{X}$	X (known from absorbance)

# Data Analysis - Equilibrium Constant

$$K = \frac{[FeSCN^{2+}]_{equilibrium}}{[Fe^{3+}]_{equilibrium} [SCN^{-}]_{equilibrium}}$$

$$K = \frac{x}{([Fe^{3+}]_{initial} - x)([SCN^{-}]_{initial} - x)}$$

where  $x = [FeSCN^{2+}]_{equilibrium}$  determine by absorption measurement

Calculate K for Solutions 8 to 12 and calculate the average K.